

# Aberrant glycosylation as a biomarker in oncology

## *Spremembe v glikozilaciji kot biomarker v onkologiji*

**Neža Gregorčič<sup>1</sup>, Shruti Jain<sup>2</sup>, Kamlesh Gidwani<sup>2</sup>**

<sup>1</sup>Univerza v Ljubljani, Medicinska fakulteta, Inštitut za biokemijo in molekularno genetiko

<sup>2</sup>University of Turku, Faculty of Technology, Department of Life Technologies, Biotechnology Unit

*Avtor za korespondenco:*

**Neža Gregorčič, Master of Science, Republika Finska**

Univerza v Ljubljani, Medicinska fakulteta, Inštitut za biokemijo in molekularno genetiko, Vrazov trg 2, 1000 Ljubljana

e-pošta: neza.gregorcic@mf.uni-lj.si

### ABSTRACT

Cancer represents a significant global health challenge, with nearly 20 million new cases and over 9 million deaths annually. Early diagnosis is a key factor in improving survival rates, yet current diagnostic methods lack sensitivity, specificity, and accessibility. Glycosylation patterns of biomarkers offer a promising target for early detection of cancer and the development of personalized treatment strategies. The europium-doped nanoparticle-based glycovariant assay represents a novel approach that leverages time-resolved fluorescence for enhanced sensitivity and glycan moiety recognition for enhanced specificity. This review explores glycosylation changes in cancer, emphasizing their potential as diagnostic tools, with a particular focus on the europium-doped nanoparticle-based glycovariant assay. The assay has demonstrated potential for identifying specific cancer-associated glycovariants in ovarian, prostate, breast, colorectal, and bladder cancers. The assay has the capacity to distinguish cancer-specific glycovariants from nonmalignant sources. However, further optimization and bigger sample cohorts are necessary to enhance the diagnostic utility of the assay.

**Key words:** glycosylation, cancer, biomarker, time-resolved fluorescence, lectins

### POVZETEK

Rak je s skoraj 20 milijoni novih primerov in 9 milijoni smrti letno znatni izziv globalnemu zdravju. Za izboljšanje stopnje preživetja je ključna zgodnja diagnoza, vendar so trenutne diagnostične metode, kar se tiče občutljivosti, specifičnosti in dostopnosti, pomanjkljive. Spremembe v glikozilaciji biomarkerjev nudijo obetavno možnost za zgodnje odkritje in izdelavo personaliziranih strategij zdravljenja. Analiza glikovariant s pomočjo z evropijem polnjenih nanodelcev je nov pristop, ki uporablja fluorescenco s časovnim zamikom za povečano občutljivost in prepoznavanje sladkorne komponente molekul za večjo specifičnost. Članek se osredotoča na spremembe glikozilacij pri raku s poudarkom na njihovem diagnostičnem potencialu in metodah, v ospredju je s pomočjo z evropijem polnjenih nanodelcev. Ta metoda kaže visok potencial za identifikacijo določenih glikovariant, ki jih povezujeemo z rakom jajčnikov, prostate, dojke, črevesja in mehurja. Zmožnost razlikovanja za raka specifičnih glikovariant in glikovariant nemalignega izvora odpira vrata za izboljšano obravnavo onkoloških pacientov. Za večjo uporabnost metode pa so še vedno potrebne optimizacije in analiza z večjim številom vzorcev.

**Ključne besede:** glikozilacija, rak, biomarker, fluorescenca s časovnim zamikom, lektini



## INTRODUCTION

With nearly 20 million new cases and over 9 million deaths reported annually, cancer remains a pressing global problem (1). The prognosis for most cancers could be markedly enhanced through the implementation of early disease intervention, which hinges on the ability to diagnose the disease at an early stage. However, current diagnostic and monitoring strategies often fall short in terms of sensitivity, specificity, and accessibility. While traditional methods, such as imaging techniques and tissue biopsies, have proven invaluable, there are limitations to their effectiveness. They are invasive, costly, and susceptible to sampling errors. Moreover, existing biomarkers, though informative in certain contexts, may lack the necessary specificity or sensitivity for early detection or accurate monitoring of disease progression (2–4).

A biomarker can be defined as “a characteristic that is objectively measured and evaluated as an indicator of normal biological processes, pathogenic processes, or pharmacologic responses to a therapeutic intervention” (5). A cancer biomarker is considered clinically useful if it is significantly elevated compared to levels in age-matched healthy and benign controls and/or reflective of the disease progression/regression, and is present in easily obtained bodily fluids, such as serum, plasma, and urine (6, 7). The characteristics of an ideal biomarker are largely dependent on its intended use. For screening purposes, the biomarker should be highly specific, reference values should be well known, it should be cost-effective, and there should be implications for therapeutic or lifestyle interventions. Conversely, a diagnostic biomarker should be highly sensitive, cost-effective, have therapeutic implications, and elucidate the molecular pathogenesis of the disease. Finally, an ideal biomarker for prognosis and treatment monitoring should be highly specific, with low intra-individual variability, it should respond to therapy, be predictive of treatment effectiveness, and add to already known prognostic indicators or models (8). The majority of cancer biomarkers currently used in clinical practice lack the necessary performance characteristics to be of diagnostic value. Consequently, the diagnostic use of biomarkers depends largely on emerging clinical symptoms, family history, or genetic susceptibility. The development of cancerspecific biomarkers that are easily accessible and sensitive is required to improve screening, early diagnosis, monitoring, and prediction of disease progression and treatment (6, 7, 9).

In recent years, there has been a growing interest in utilizing cancer-specific glycosylation patterns as diagnostic and prognostic biomarkers. These distinct glycan signatures have the potential to facilitate early detection, enable the development of personalized treatment strategies, and facilitate the assessment of disease progression (10). This article provides an overview of glycosylation changes in cancer, emphasizing their potential as diagnostic tools, with a particular focus on the europium-doped nanoparticle-based glycovariant (EuNP-based GV) assays for the detection of glycosylated cancer biomarkers.

## GLYCOSYLATION

Glycans are saccharide chains, which can either be free or in the form of glycoconjugates. As glycoconjugates, they are attached to lipids or proteins. This attachment of sugars to proteins or lipids occurs during a posttranslational process known as glycosylation, which occurs in the endoplasmic reticulum and Golgi apparatus (11). The process is done in several steps, which depend on substrate and enzyme availability (11, 12). Glycosylation begins with the determination of proteins to be glycosylated and the transfer of the monosaccharide to an amino acid. This is a part of the initiation step, which occurs in the endoplasmic reticulum. The exceptions to this are two Oglycosylations, N-acetylgalactosamine-type and Xylitol-type, which are initiated in the Golgi apparatus, and O-linked-N-acetylglucosaminylation, which occurs in the cytosol and nucleus. The initiation step is followed by glycotransferases further adding monosaccharides to growing oligosaccharide chains, resulting in core extension, elongation, and branching, as well as the final capping of glycans (12).

Most of the membrane and secreted proteins are glycosylated. Their variety and abundance directly reflect their importance and many functions. Their roles can be divided into three general categories: structural and modulatory properties, specific recognition by other molecules, and the molecular mimicry of host glycans. The first category consists of structural and modulatory properties, which include protective, stabilizing, organizational, and barrier functions. For example, glycocalyx acts as a physical barrier, shielding cellular contents from proteolytic recognition, antibody binding, and microbial attachment. Glycosylation also modulates protein interactions and protects

»


important molecules by binding them to glycosaminoglycan chains within the extracellular matrix, preventing the diffusion of factors away from the site and prolonging their active lives. The second category is specific recognition, which involves intrinsic interactions between glycans and glycan-binding proteins on cell surfaces. These interactions are crucial for cell-cell and cell-matrix interactions. An example is selectins, which recognize glycans on ligands to mediate the binding of blood and vascular cells in various situations. Finally, the third category, molecular mimicry of host glycans, refers to glycans that are specifically bound by various viruses, bacteria, and parasites, and glycans, targeted by toxins. These glycans play numerous roles. For example, some organisms can mask or modify glycans typically recognized by microorganisms or toxins, thus preventing their binding or damage (13).

### Glycosylation changes in cancer

In the context of tumor progression, several glycan aberrations have been identified as pivotal pathophysiological events. These include the loss or overproduction of some glycans, the increased production of truncated glycans, and the synthesis of novel glycans (14–16). Such alterations stem from two main mechanisms: incomplete synthesis and neo-synthesis. Incomplete synthesis often occurs in early cancer stages, where normal synthesis of complex

glycans is impaired, resulting in the biosynthesis of truncated structures. On the other hand, neo-synthesis is more prevalent in advanced cancers, where cancer cells acquire genes involved in glycan biosynthesis, leading to abnormal glycosylation patterns (10). These changes in the glycosylation pathway may occur in the expression of specific glycotransferases, glycosidases, molecular chaperones, or lysosomal hydrolases. Another possibility is the mislocalization of the glycosyltransferases or a change in their nature. Glycosylation changes may also occur due to altered availability of protein substrates or nucleotide sugars and 3-phosphoadenosine-5-phosphosulfate, activity of nucleotide sugar and monosaccharide transporters, glycoprotein turnover kinetics, and pH of endoplasmic reticulum and Golgi apparatus (17).

Glycosylation aberrations that frequently occur in cancer include an increase in sialyl-Lewis structures, abnormal core fucosylation, increased sialylation, elevated N-glycan branching or exposure of the mucin-type O-glycan, and formation of Thomsen-Friedenreich (T) and Thomsennoveau (Tn) antigens. These changes play a role in many oncogenic processes, such as evasion of cell division checkpoints, evasion of death signals and immune surveillance, and migration to metastatic sites (10, 16, 18). Cancer-associated glycans are shown in Figure 1. >>



**Figure 1: Cancer-associated glycans.** The figure shows glycans commonly found in cancers. N-acetylglucosamine (GlcNAc), Mannose (Man), Galactose (Gal), N-Acetylneuraminic Acid (Neu5Ac), Asparagine (Asn), Serine (Ser), Threonine (Thr). Created with BioRender.com.

**Slika 1: Z rakom povezani glikani.** Na sliki so prikazani glikani, pogosto prisotni pri rakavih obolenjih. N-acetilglukozamin (GlcNAc), manoza (Man), galaktoza (Gal), N-acetilneuraminska kislina (Neu5Ac), asparagin (Asn), serin (Ser), treonin (Thr). Slika je bila pripravljena s spletnim orodjem BioRender.com.



## Glycosylation of extracellular vesicles

Extracellular vesicles (EVs), membranous particles with a size range of 40 to 1000 nm, enclosed by a phospholipid bilayer, are secreted by various cells and represent a potential analyte for diagnosis and treatment monitoring of cancer (19, 20). EVs carry clinically relevant information, as they reflect the cell condition from which they originate. Furthermore, their production is increased in proportion to tumor size and growth rate (7). Moreover, since they can be found in bodily fluids like urine, plasma, and saliva, they can be obtained non-invasively (19). Similarly to free cancer-specific glycoproteins, EVs also present a unique set of glycan epitopes on their membrane (7).

Cancer-derived EVs have been found to display a tumor-associated glycocalyx, which could serve as crucial markers for the detection and isolation of EVs, and notably, for malignant and nonmalignant distinction. These changes have been observed in the early stages of the disease, indicating a potential use for early diagnosis (7, 21). Some of the reported glycocalyx aberrations include high-mannose and complex type N-glycans, polylactosamine, and sialylated glycans, specifically  $\alpha$ 2,6linked sialic acids (21).

## CANCER-SPECIFIC GLYCOSYLATION BIOMARKERS

Glycosylation changes in cancer have long been considered a hallmark of cancer and most of the currently recognized and used cancer biomarkers are glycoproteins. However, they are normally monitored for their total protein level. Currently, only two glycosylated biomarkers, alpha fetoprotein-L3 and cancer antigen 19-9 (CA19-9), are defined by their glycan moiety in clinical practice. The discovery and recognition of cancer-specific GVs could improve the specificity of these cancer biomarkers (2, 10, 16, 22–24). The most commonly employed glycoproteins used as cancer biomarkers are listed in Table 1, along with reported glycan alterations and the potential utility of glycan-specific detection. Some promising glycan targets, which are not yet in clinical use, are also listed.

**Table 1. Glycoproteins as cancer biomarkers.** Some of the most used glycoprotein cancer biomarkers and some potential new cancer biomarkers, glycosylation alterations reported in cancer, and the potential of glycovariant detection.

**Tabela 1. Glikoproteini kot onkološki biomarkerji.** Najpogostejši glikoproteini, ki se uporabljajo kot onkološki biomarkerji v praksi, in nekateri potencialni novi biomarkerji, za raka specifične spremembe v glikolizaciji in potencialna uporaba glikovariant.

Biomarker	Cancer type	Clinical use	Glycosylation alterations found in cancer	Potential utility of GV detection	Ref.
AFP-L3	Hepatic	Diagnosis, prognosis	Core-fucosylation	Already in clinical use	(27)
CA19-9	Pancreatic	Diagnosis, prognosis	Sialyl-Lewis A structure	Already in clinical use	(27)
CA15-3	Breast	Monitoring of disease progression and response to therapy	Core 1-based O-glycosylation in place of the normal core 2based O-glycosylation, Tn and T-sialylation, and 3O sulphated or 3-sialylated core 1 and extended core 1	Increased specificity and sensitivity	(68)
CA125	Ovarian	Diagnosis, monitoring treatment response, detecting recurrence	Increased core-fucosylated biantennary mono-sialylated glycans, high-mannose type, and complex type N-glycan. Decreased bisecting biantennary and nonfucosylated glycans. Core 1 and 2 type glycans with branched core 1 antennae	Increased sensitivity and specificity	(56, 57)



Biomarker	Cancer type	Clinical use	Glycosylation alterations found in cancer	Potential utility of GV detection	Ref.
CEA	Colorectal	Monitoring treatment response, detecting recurrence	Branched N-glycans, higher density of O-glycans, incomplete glycan synthesis, increased sialylation and fucosylation	Screening	(69)
PSA	Prostate	Diagnosis, monitoring treatment response, detecting recurrence	Changes in core fucosylation and sialylation, highly branched N-glycans, low molecular weight PSA without N-glycans, and Lewis X antigen	Improved specificity	(66, 71)
Hp	Hepatic	Prognosis	Increased bi-fucosylation	Diagnosis	(27, 72)
	Ovarian	Prognosis	Increased fucosylation	Diagnosis	
Fetuin-A	Hepatic		Increased fucosylation	Diagnosis, prognosis	(27)
	Cholangio				
Total serum/plasma profiles	Hepatic		Increased multi-antennary glycans with fucose residues	Diagnosis	
	Gastric		Decreased core-fucosylated glycans; increased hybrid and multi-branched structures and decreased monoantennary, galactose, bisecting type, and core fucose structures	Diagnosis, monitoring progression	
	Pancreatic		Increased core-fucose residues	Diagnosis	
	Breast		Increased sialylation, branching, and outer-arm fucosylation and decreased high-mannosylated and biantennary core-fucosylated glycans	Diagnosis	
IgG	Hepatic		Decreased galactosylation; increased core-fucosylation	Diagnosis, prognosis	
	Colorectal		Decreased galactosylation; decreased galactosylation and sialylation; increased core-fucosylation of neutral glycans and decreased corefucosylation of sialylated glycans	Diagnosis	
	Gastric		Decreased galactosylation	Diagnosis	
	Lung		Decreased galactosylation	Diagnosis	
	Ovarian		Decreased galactosylation	Diagnosis	
AIAT	Lung		Increased galactosylation, fucosylation, and poly-LacNAc structures	Diagnosis	
	Hepatic		Increased fucosylation	Diagnosis	




Glycovariant (GV); References (Ref.); Alpha-fetoprotein (AFP-L3); Cancer antigen 15-3 (CA153); Cancer antigen 19-9 (CA19-9); Cancer antigen 125 (CA125); Carcinoembryonic Antigen (CEA); Prostate-Specific Antigen (PSA); Haptoglobin (Hp); Immunoglobulin G (IgG);  $\alpha_1$ Antitrypsin (AIAT);  $\alpha_1$ -Acid glycoprotein (AGP); Ceruloplasmin (CP)

## Europium-doped nanoparticles-based glycovariant assay

Several advanced techniques are used to analyze glycoproteins to target glycan alterations in cancer, each with distinct advantages and limitations. Detailed glycoprotein characterization methods, such as high-performance liquid chromatography, mass spectrometry, and nuclear magnetic resonance spectroscopy, provide high sensitivity and specificity for detailed structural and quantitative analysis but require complex preparation, skilled personnel, expensive equipment, and a lot of time. Thus these methods, while useful for novel biomarker discovery, are not suitable for routine clinical use (25–27). To this end, a range of lectin-based assays has been developed, including microarrays, enzyme-linked lectin assays, lectin immunosorbent assays, and affinity chromatography. Their key advantage is the ability to differentiate between vario-

us glycan structures without the removal of the sugar component from glycoproteins or the addition of fluorescent labels. Although these methods show promise for the detection of cancer glycoproteins because of their simplicity and suitability for high-throughput screening, they exhibit several disadvantages, such as the need to deglycosylate antibodies used for capture or detection, and lower sensitivity and specificity due to the binding of lectins to glycans expressed by unrelated proteins (25, 27).

The EuNP-based GV assay is a sandwich-like immunoassay. In this assay, immobilized capture antibodies recognize and bind the protein moiety of the analyte, while lectins or glycan-specific antibodies coated on EuNPs recognize and bind the cancer-specific glycan moiety. The amount of glycoprotein is then measured as time-resolved fluorescence (28). This principle is shown in Figure 2.



**Figure 2: Europium-doped nanoparticle-based glycovariant assay.** In this assay, immobilized capture antibodies recognize and bind the protein moiety of the analyte, while lectins or glycan-specific antibodies coated on nanoparticles recognize and bind the cancer-specific glycan moiety. The amount of glycoprotein is then measured as time-resolved fluorescence. Created with BioRender.com.

**Slika 2: Analiza glikovariant s pomočjo z evropijem polnjenih nanodelcev.** Pri tej metodi imobilizirana protitelesa prepoznajo in vežejo beljakovinski del analita, medtem ko lektini ali za glikane specifična protitelesa, vezana na zunanosti nanodelcev, prepoznajo in vežejo za raka specifični glikanski del. Količina glikoproteina se nato izmeri kot fluorescenca s časovnim zamikom. : Slika je bila pripravljena s spletnim orodjem BioRender.com.





The use of different capture and tracer (captured glycoprotein<sup>tracer</sup>) combinations in EuNP-based GV assays has shown promise for the specific detection of cancer-associated glycovariants across multiple cancer types, including prostate, breast, colorectal, and bladder (28–38). Furthermore, this method has shown potential for the detection of glycosylation changes on the surface of EVs. It has been proposed that cancer-specific EVs could be effectively characterized using capture antibodies and EuNP-linked lectins and antibodies targeting specific proteins or glycan structures known to be present on the EV surface, such as tetraspanins (e.g., CD9, CD63, CD81) and integrins (20,39–42).

## COMPONENTS OF THE EUROPIUM-DOPED NANOPARTICLE-BASED GLYCOVARIANT ASSAY

The performance of the EuNP-based glycovariant assay hinges on several essential equipment components and materials. Firstly, the assay requires the use of fluorescent NPs, specifically (Eu<sup>3+</sup>)-chelate-doped NPs, which exhibit time-resolved fluorescence. Consequently, a time-resolved fluorescence microplate reader is required to measure the signal produced. Secondly, lectins or glycan-specific antibodies must be coated on the aforementioned NPs so that they can bind to the glycan moiety of the target glycoprotein. Thirdly, biotinylated capture antibodies must be immobilized on a streptavidin-coated low-fluorescence microtitration plate in order to bind the protein moiety of the target glycoprotein. Additionally, some basic materials, such as wash and assay buffers, are required to maintain proper conditions during the assay. Lastly, conventional laboratory equipment, such as plate shakers for incubation and centrifuges for sample preparation, is needed (28–30, 32, 35, 43).

### Lectins

Lectins are proteins of nonimmune origin that specifically recognize and reversibly bind glycans without changing their structure (15, 44). They can form bonds with

monosaccharides, oligosaccharides, and sugar residues in complex molecules. The carbohydrate-lectin interaction depends on many factors, such as lectin valence, the structure of binding sites, and their spatial arrangement (44).

As an analytical tool, the most important aspect of lectins is their specificity. This characteristic is determined by carbohydrate residues in glycans (44, 45). Each lectin exhibits a specific affinity for a particular carbohydrate structure (45). Carbohydrate residues may be located either at the terminus of a carbohydrate chain or within the chain itself. Some lectins can recognize and bind terminal residues, while others can bind those within a chain. Lectins are often specific to the definite anomeric form of carbohydrate molecules ( $\alpha$ - or  $\beta$ -anomer) and they differentiate between sequences of carbohydrate residues (44).

Another crucial aspect of lectins' analytical utility is their affinity for the analyte. Unlike antibodies, which exhibit high affinity with dissociation constants in the range of  $10^{-8}$  to  $10^{-12}$  M, lectins typically have lower affinity, with dissociation constants ranging between  $10^{-4}$  and  $10^{-7}$  M. However, when lectins interact with polysaccharides through polyvalent interactions, their affinities can be comparable to those observed in immune reactions (38, 44).

### Glycan-binding antibodies

Antibodies possess diverse structures, enabling them to bind specifically and with high affinity to a wide range of targets. However, their effectiveness in detecting carbohydrates is limited compared to proteins due to various factors, such as the less T-cell dependent nature of anti-glycan immune responses and the preference for larger, non-carbohydrate portions of biomolecules. Despite these challenges, efforts have been made to develop anti-glycan antibodies against various targets including cancer epitopes, such as the monoclonal anti-Sialyl-Tn (anti-STn) antibody. However, the process of generating these antibodies is often tedious due to the lack of accessibility of pure glycans. The antibodies themselves can also suffer from issues such as specificity and affinity (46).

### Europium-doped nanoparticles

Polystyrene NPs, employed as labels in the EuNP-based GV assay, have a diameter of 107 nm and contain approximately 30 000 chelated Eu<sup>+3</sup> ions per particle. Due to the

»



high concentration of  $\text{Eu}^{+3}$  chelates in each NP, the signal is amplified by a factor of 30 000 compared to a single  $\text{Eu}^{+3}$  chelate. The coupling of several antibodies or lectins on the surface of a single NP enables the creation of a high density of binding sites on the surface of NPs, thereby increasing their avidity. The combination of high avidity and signal amplification significantly enhances the sensitivity of the assay (6, 47).

### Time-resolved fluorescence

Fluorescence is a phenomenon whereby a fluorophore, a molecule with the capacity to fluoresce, absorbs electromagnetic radiation, typically from ultraviolet or visible light, and subsequently emits a photon of lower energy, which implies a smaller frequency and a longer wavelength. The distinction between the absorbed and emitted energy is referred to as the Stokes shift (48). Fluorometric measurement presents certain challenges, including autofluorescence and nonspecific fluorescence, which can result in decreased sensitivity (49).

Lanthanides, such as Eu, and their chelates, exhibit distinctive fluorescent properties, including long fluorescence decay times (up to over 2 ms when in a complex), sharp emission peaks, and large Stokes shifts (50, 51). These properties render them an optimal label for detection with time-resolved fluorescence. The application of time-resolved fluorescence measurement has been shown to reduce the background interference of biological matrix fluorescence and to improve the sensitivity of such assays (51).

### Samples

The sensitivity of the EuNP-based GV assay allows for its performance on small sample volumes (a few microliters), which is an advantage over traditional enzyme-linked immunosorbent assays (52). The types of samples required for these assays are typically bodily fluids in which cancer markers are present. Commonly used samples include serum, plasma, and more specialized fluids such as ovarian cyst fluid, placental homogenate, and ascites fluid (28, 32, 34, 36, 52, 53). The choice of sample depends on the specific cancer being investigated (52).

## Ovarian cancer

Ovarian cancer is one of the most prevalent gynecological cancers, with over 300 000 reported new cases worldwide according to Globocan 2022 statistics (1). The primary challenge associated with this cancer is a high mortality rate, which is largely attributable to late-stage diagnosis in the majority of cases. Currently, diagnosis is based on serum cancer antigen 125 (CA125) concentration in symptomatic women. If values are elevated, an ultrasound of the abdomen and pelvis is used to confirm the diagnosis (54).

CA125 is a membrane-bound mucin – a highly O-glycosylated glycoprotein. O-glycans account for over 50% of its mass. The high amounts of glycosylation are due to tandem repeats with high amounts of serines and threonines creating many glycosylation sites, which leads to many subspecies. As a serum biomarker, it is typically employed for the diagnosis and monitoring of epithelial ovarian cancer. However, elevated concentrations have been observed in several other conditions, including endometriosis, liver disease, ovarian cysts, pregnancy, ovulatory cycles, heart failure, and some other malignant conditions (55). The biomarker's extensive glycosylation suggests the potential for increased specificity if cancer-specific glycovariants can be distinguished. Elevated concentrations of core-fucosylated biantennary monosialylated glycans in serum samples of cancer patients have been reported. In contrast, levels of mostly bisecting biantennary and non-fucosylated glycans have been found to be decreased. Furthermore, a prevalence of high-mannose type and complex type Nglycan structures has been observed. In addition, O-glycans in ovarian cancer cell lines, such as OVCAR3, are predominantly composed of core 1 and 2 type glycans with branched core 1 antennae (56, 57).

### Europium-doped nanoparticles-based glycovariant assay in ovarian cancer

The EuNP-based GV assay was used to identify lectins that can discriminate between CA125 from epithelial ovarian cancer and non-cancer sources. It was found that macrophage galactosyltransferase lectin (MGL)-coated EuNPs specifically bind to the CA125 isoform produced by cancer cells. In addition, the detection of CA125 from non-malignant sources was reduced compared to the conventional immunoassay (58). Similar observations were made with anti-STn monoclonal antibodies coated on  $\text{Eu}^{+3}$ -doped NPs, which efficiently recognized the epithelial ovarian cancer- >>

rassociated CA125 isoform produced by OVCAR-3. The assay also detected this CA125 GV in serum with a sensitivity of 73.3%, while 95% of endometriosis cases remained undetected (59). The performance of the CA125<sup>STn</sup> assay is particularly impressive at high specificity (90% or greater), with improvements of 10.6% and 8.9% over conventional CA125 immunoassay with and without borderline cases. An even more dramatic improvement was seen in postmenopausal cases, with 17.2% higher detection at 90% specificity compared to the conventional immunoassay (60).

Both MGL and anti-STn coated NPs were also evaluated as a potential tool in the differential diagnosis of pelvic masses in a larger sample cohort; 549 patients with epithelial ovarian cancer, benign ovarian tumors, and endometriosis. Both GV assays demonstrated superior area under the ROC curve (AUC) and lower false positive rates compared to the conventional CA125 immunoassay. The sensitivity of the CA125<sup>STn</sup> assay was significantly better than the conventional immunoassay, with 85% vs. 74% sensitivity at 90% specificity (61). Both GVs were also studied to determine their potential for monitoring treatment and follow-up in high-grade serous ovarian cancer. The study showed promising results for the CA125<sup>STn</sup> assay in differentiating between low and high tumor burden. Patients with complete cytoreduction were found to have significantly lower STn glycoprotein levels at diagnosis than those with suboptimal cytoreduction. In contrast, patients with suboptimal cytoreduction could not be differentiated using the conventional CA125 immunoassay. The nadir value of CA125<sup>STn</sup> was predictive of the patients' progression-free survival. The CA125<sup>STn</sup> test improved the detection of disease recurrence, showing a higher biomarker increase in 80.0% of patients and an earlier increase in 37.0% of patients (62).

While CA125<sup>STn</sup> appears to be the most promising combination to date, other EuNP-based GV assays have also shown some potential in discriminating epithelial ovarian cancer, such as CA153<sup>STn</sup>, CA15-3<sup>Tn</sup>, integrin  $\alpha 3^{\text{STn}}$ , and CD63<sup>STn</sup> (39, 60).

## PROSTATE CANCER

With about 1.5 million newly reported cases in 2022, prostate cancer ranks as the most common cancer in males world-

wide (1). The most commonly used biomarker and one of the oldest biomarkers assisting in prostate cancer diagnosis is the serum prostate-specific antigen (PSA). PSA is a kallikrein glycoprotein with a single N-oligosaccharide chain attached to Asn-54 produced by prostate epithelial cells (63, 64). An elevation in serum PSA levels is typically indicative of disruption to the prostate basement membrane, which is normally caused by cancer cells. However, it can also be caused by benign prostatic hyperplasia (BPH), prostatitis, or manipulation of the prostate. A lack of specificity frequently results in the overdiagnosis of patients and the administration of unnecessary treatments (63, 65). A significant proportion of PSA's mass is constituted by heavy O-glycosylation. Studies have reported alterations in glycosylation patterns, including increased core fucosylation and sialylation, in prostate cancer. Detection of these aberrant glycovariants could improve discrimination between malignant and benign conditions (66).

### Europium-doped nanoparticles-based glycovariant assay in prostate cancer

The EuNP-based GV assay was employed to specifically identify fucosylated PSA. The glycan moiety was bound by Aleuria aurantia lectin as a proof-of-concept study on a small cohort of 36 tissue samples (11 benign and 25 cancerous) and 42 urine samples (16 cancerous, 15 BPH, and 11 young males). The findings revealed a markedly elevated PSA fucosylation level in malignant tissue in comparison to benign tissue. In urine samples, the EuNP-based GV assay demonstrated superior discrimination between the three groups compared to total PSA and free PSA concentrations in urine. However, the promising results observed in urine samples did not translate into serum and plasma samples due to the high variability of the LN-CaP cell line PSA standard recovery and high background signal. EVs may also represent a promising target for the GV assay. Higher integrin  $\alpha 3$  levels were observed in more aggressive prostate cancer-derived cell lines, whereas less aggressive cell lines exhibited higher epithelial cell adhesion molecule levels (40). A further investigation into the potential of lectins for cancer biomarker detection has been conducted, with five lectins identified as having the greatest promise for detection in serum and plasma. These include mannose-binding lectin, *Trichosanthes japonica* agglutinin II, Wheat Germ Agglutinin (WGA), *Wisteria floribunda* lectin and MGL (6).

»

## BREAST CANCER

Breast cancer is the most common cancer in women, with over 2 million reported new cases in 2022 (1). Cancer antigen 15-3 (CA15-3) is normally expressed in secretory epithelial cells at the apical plasma membrane. Elevated CA15-3 is observed in the majority of breast cancer patients. Elevation of the mucin also occurs in different types of advanced adenocarcinoma, such as pancreatic, gastric, and lung adenocarcinoma (67). Alterations in the glycosylation pattern of CA15-3 have been observed in breast cancer, including the presence of core 1-based Oglycosylation in place of the normal core 2-based O-glycosylation, Tn and T-sialylation, and 3O-sulfated or 3-sialylated core 1 and extended core 1 (68).

### Europium-doped nanoparticles-based glycovariant assay in breast cancer

CA15-3 was glycoprofiled with 28 lectins to identify those that recognize the carbohydrate structural changes associated with breast cancer. The most promising tracers were MGL and WGA-doped NPs, with clinical sensitivities of 67.9% and 81.1%, respectively, and specificities of 90%. In comparison, the conventional CA15-3 assay demonstrated a sensitivity of 66.0% at 90% specificity (68). In a subsequent study, the CA15-3<sup>WGA</sup> assay was again identified as a promising assay for breast cancer detection, in conjunction with the CA125<sup>WGA</sup> assay and the CD63 immunoassay. Significant differences were observed between the case and control groups for the CA125<sup>WGA</sup>, CA15-3<sup>WGA</sup>, and CD63 immunoassay measurements. Furthermore, significant improvements were reported in the AUCs of the new GV assays, particularly in combination with CD63 immunoassay, as compared to the conventional immunoassay (41).

## COLORECTAL CANCER

Colorectal cancer is the third most common cancer worldwide, with over a million reported new cases in 2022 (1). Screening programs that identify precancerous polyps or diagnose colorectal cancer using colonoscopy tend to be successful for the 80% of colorectal cancers that develop from these polyps (69). Biomarkers can prove particularly

useful in monitoring colorectal cancer treatment. Currently, the most widely used biomarker for this purpose is carcinoembryonic antigen (CEA). However, high concentrations of CEA also appear in benign conditions such as hepatitis, pancreatitis, inflammatory bowel disease, and Crohn's disease. Other biomarkers can assist with diagnosis, postoperative follow-up, and treatment monitoring. These include CA19-9, CA125, cancer antigen 72-4 (CA72-4), and tumor-associated glycoprotein 72 (TAG 72) (40, 50). Several glycosylation processes occur during malignant transformation in colorectal cancer, including increased branching of N-glycans, higher density of O-glycans, incomplete glycan synthesis, glycan neo-synthesis, increased sialylation, and increased fucosylation (69).

### Europium-doped nanoparticles-based glycovariant assay in colorectal cancer

A small cohort of 22 colorectal cancer patients was evaluated, comprising 11 early-stage (I–II) and 11 late-stage (III–IV) cases. CA125<sup>WGA</sup>, CA19-9<sup>Ma696</sup>, and CA19-9<sup>Con-A</sup> assays demonstrated high sensitivity and specificity in detecting late-stage samples, whereas the CD151<sup>CD63</sup> assay exhibited limited performance in recognizing early-stage colorectal cancer. CA125<sup>WGA</sup> identified 11 out of 22 colorectal cancers, CA19-9<sup>Ma696</sup> recognized 9 out of 22 colorectal cancers, CA19-9<sup>Con-A</sup> detected 11 out of 22 colorectal cancers, and CD151<sup>CD63</sup> assay detected 8 out of 11 early-stage colorectal samples. A combination of all three assays may be useful for the post-treatment follow-up, for the detection of early relapse, and for the screening of colorectal cancer. However, a larger sample cohort is necessary for the evaluation of the assays (42).

## BLADDER CANCER

Bladder cancer is the second most common urological cancer, with over 6 hundred thousand new cases reported in 2022 (1). The current diagnostic approach for bladder cancer is based on cystoscopy and urine cytology. While cystoscopy is an invasive and expensive procedure, cytology lacks sensitivity to low-grade tumors, necessitating the development of a targeted, less invasive, more sensitive, and highly specific approach for bladder cancer diagnosis (70).

## Europium-doped nanoparticles-based glycovariant assay in bladder cancer

One potential avenue for improving the diagnosis process is the utilization of Ulex europaeus agglutinin-I (UEA-I)-coated NPs in a GV integrin  $\alpha 3^{UEA-I}$  assay. The assay was conducted on a limited cohort of urine samples, comprising 13 patients with bladder cancer and 9 patients with BPH. The assay demonstrated encouraging outcomes, with a statistically significant differentiation between cancer and benign samples. However, further research is necessary (20).

## CONCLUSION

In conclusion, cancer remains a major global health problem and early diagnosis is crucial for improving survival rates. Current diagnostic methods have limitations, and there is a need for easily accessible and sensitive cancer biomarkers. Glycosylation patterns have shown promise as diagnostic and prognostic biomarkers, and the EuNP-based GV assay has demonstrated potential in detecting glycosylated cancer biomarkers. However, further research and development are needed to validate these biomarkers and improve their sensitivity and specificity.

## REFERENCES

1. Ferlay J, Ervik M, Lam F, Laversanne M, Colombet M, Mery L, et al. Global Cancer Observatory: Cancer Today [Internet]. Lyon, France: International Agency for Research on Cancer. [accessed on 20. 5. 2024]. Available from: <https://gco.iarc.fr/today/en>
2. He K, Baniasad M, Kwon H, Caval T, Xu G, Lebrilla C, et al. Decoding the glycoproteome: a new frontier for biomarker discovery in cancer. *J Hematol Oncol.* 2024;17(1).
3. Crosby D, Bhatia S, Brindle KM, Coussens LM, Dive C, Emberton M, et al. Early detection of cancer. *Science.* 2022;375(6586).
4. Martins I, Ribeiro IP, Jorge J, Gonçalves AC, Sarmiento-Ribeiro AB, Melo JB, et al. Liquid biopsies: Applications for cancer diagnosis and monitoring. *Genes (Basel).* 2021;12(3):1–20.
5. Atkinson AJ, Colburn WA, DeGruttola VG, DeMets DL, Downing GJ, Hoth DF, et al. Biomarkers and surrogate endpoints: Preferred definitions and conceptual framework. *Clin Pharmacol Ther.* 2001;69(3):89–95.
6. Gidwani K, Kekki H, Terävä J, Soukka T, Sundfeldt K, Pettersson K. Nanoparticle-aided glycovariant assays to bridge biomarker performance and ctDNA results. *Mol Aspects Med.* 2020;72.
7. Grzesik K, Janik M, Hoja-Lukowicz D. The hidden potential of glycomarkers: Glycosylation studies in the service of cancer diagnosis and treatment. *Biochim Biophys Acta Rev Cancer.* 2023;1878(3).
8. Burlina A, Brand E, Hughes D, Kantola I, Krämer J, Nowak A, et al. An expert consensus on the recommendations for the use of biomarkers in Fabry disease. *Mol Genet Metab.* 2023;139(2).
9. Sarhadi VK, Armengol G. *Molecular Biomarkers in Cancer.* Biomolecules. 2022;12(8).
10. Pinho SS, Reis CA. Glycosylation in cancer: Mechanisms and clinical implications. *Nat Rev Cancer.* 2015;15(9):540–55.
11. Reily C, Stewart TJ, Renfrow MB, Novak J. Glycosylation in health and disease. *Nat Rev Nephrol.* 2019;15(6):346–66.
12. Schjoldager KT, Narimatsu Y, Joshi HJ, Clausen H. Global view of human protein glycosylation pathways and functions. *Nat Rev Mol Cell Biol.* 2020;21(12):729–49.
13. Varki A, Gagneux P. Biological Functions of Glycans. In: Varki A, Cummings R, Esko J, editors. *Essentials of Glycobiology.* 3rd ed. New York: Cold Spring Harbor Laboratory Press; 2017.
14. Bellis S, Reis C, Varki A, Kannagi R, Stanley P. Glycosylation Changes in Cancer. In: Varki A, Cummings R, Esko J, et al., editors. *Essentials of Glycobiology.* 4th ed. New York: Cold Spring Harbor Laboratory Press; 2022.
15. Fuster MM, Esko JD. The sweet and sour of cancer: Glycans as novel therapeutic targets. *Nat Rev Cancer.* 2005;5(7):526–42
16. Hakomori S. Tumor Malignancy Defined by Aberrant Glycosylation and Sphingo(glyco)lipid Metabolism. *Cancer Res.* 1996;56(23):5309–5318
17. Stowell SR, Ju T, Cummings RD. Protein glycosylation in cancer. *Annu Rev Pathol.* 2015;10:473–510
18. Dimitroff CJ. I-branched carbohydrates as emerging effectors of malignant progression. *Proc Natl Acad Sci U S A.* 2019;116(28):13729–37
19. Moghassemi S, Dadashzadeh A, Sousa MJ, Vlieghe H, Yang J, León-Félix CM, et al. Extracellular vesicles in nanomedicine and regenerative medicine: A review over the last decade. *Bioact Mater.* 2024;36:126–56
20. Islam MK, Syed P, Dhondt B, Gidwani K, Pettersson K, Lamminmäki U, et al. Detection of bladder cancer with aberrantly fucosylated ITGA3. *Anal Biochem.* 2021;628.
21. Martins AM, Ramos CC, Freitas D, Reis CA. Glycosylation of cancer extracellular vesicles: Capture strategies, functional roles and potential clinical applications. *Cells.* 2021;10(1):1–26
22. Meany DL, Chan DW. Aberrant glycosylation associated with enzymes as cancer biomarkers. *Clin Proteomics.* 2011;8(1)
23. Hashim OH, Jayapalan JJ, Lee CS. Lectins: An effective tool for screening of potential cancer biomarkers. *PeerJ.* 2017;5:e3784
24. Das S, Dey MK, Devireddy R, Gartia MR. Biomarkers in Cancer Detection, Diagnosis, and Prognosis. *Sensors.* 2024;24(1).
25. Kekki H, Peltola M, van Vliet S, Bangma C, van Kooyk Y, Pettersson K. Improved cancer specificity in PSA assay using Aleuria aurantia lectin coated Eu-nanoparticles for detection. *Clin Biochem.* 2017;50(1–2):54–61.
26. Haga Y, Ueda K. Glycosylation in cancer: its application as a biomarker and recent advances of analytical techniques. *Glycoconj J.* 2022;39(2):303–13.
27. Wang M, Zhu J, Lubman DM, Gao C. Aberrant glycosylation and cancer biomarker discovery: A promising and thorny journey. *Clin Chem Lab Med.* 2019;57(4):407–16.





28. Gidwani K, Nadeem N, Huhtinen K, Kekki H, Heinosaalo T, Hynninen J, et al. Europium Nanoparticle-Based Sialyl-Tn Monoclonal Antibody Discriminates Epithelial Ovarian Cancer-Associated CA125 from Benign Sources. *J Appl Lab Med*. 2019;4(3):299–310.
29. Jain S, Nadeem N, Ulfenborg B, Mäkelä M, Ruma SA, Terävä J, et al. Diagnostic potential of nanoparticle aided assays for MUC16 and MUC1 glycovariants in ovarian cancer. *Int J Cancer*. 2022;151(7):1175–84.
30. Jain S, Parimelazhagan Santhi P, Vinod R, Afrin Ruma S, Huhtinen K, Pettersson K, et al. Aberrant glycosylation of  $\alpha 3$  integrins as diagnostic markers in epithelial ovarian cancer. *Clin Chim Acta*. 2023;543.
31. Salminen L, Nadeem N, Rolfsen AL, Dørum A, Laajala TD, Grønman S, et al. Exploratory Analysis of CA125-MGL and-STn Glycoforms in the Differential Diagnostics of Pelvic Masses. *J Appl Lab Med*. 2020;5(2):263–72.
32. Gidwani K, Huhtinen K, Kekki H, Van Vliet S, Hynninen J, Koivuviita N, et al. A nanoparticle-lectin immunoassay improves discrimination of serum CA125 from malignant and benign sources. *Clin Chem*. 2016;62(10):1390–400.
33. Terävä J, Verhassel A, Botti O, Islam MK, Leivo J, Wittfooth S, et al. Primary breast cancer biomarkers based on glycosylation and extracellular vesicles detected from human serum. *Cancer Rep*. 2022;5(8).
34. Islam MK, Syed P, Dhondt B, Gidwani K, Pettersson K, Lamminmäki U, et al. Detection of bladder cancer with aberrantly fucosylated ITGA3. *Anal Biochem*. 2021;628.
35. Salminen L, Nadeem N, Jain S, Grønman S, Carpén O, Hietanen S, et al. A longitudinal analysis of CA125 glycoforms in the monitoring and follow up of high grade serous ovarian cancer. *Gynecol Oncol*. 2020;156(3):689–94.
36. Vinod R, Mahran R, Routila E, Leivo J, Pettersson K, Gidwani K. Nanoparticle-aided detection of colorectal cancer-associated glycoconjugates of extracellular vesicles in human serum. *Int J Mol Sci*. 2021;22(19).
37. Bayoumy S, Hyytiä H, Leivo J, Talha SM, Huhtinen K, Poutanen M, et al. Glycovariant-based lateral flow immunoassay to detect ovarian cancer-associated serum CA125. *Commun Biol*. 2020;3(1).
38. Syed P, Gidwani K, Kekki H, Leivo J, Pettersson K, Lamminmäki U. Role of lectin microarrays in cancer diagnosis. *Proteomics*. 2016;16(8):1257–65.
39. Jain S, Parimelazhagan Santhi P, Vinod R, Afrin Ruma S, Huhtinen K, Pettersson K, et al. Aberrant glycosylation of  $\alpha 3$  integrins as diagnostic markers in epithelial ovarian cancer. *Clin Chim Acta*. 2023;543.
40. Islam MK, Syed P, Lehtinen L, Leivo J, Gidwani K, Wittfooth S, et al. A Nanoparticle-Based Approach for the Detection of Extracellular Vesicles. *Sci Rep*. 2019;9(1).
41. Terävä J, Verhassel A, Botti O, Islam MK, Leivo J, Wittfooth S, et al. Primary breast cancer biomarkers based on glycosylation and extracellular vesicles detected from human serum. *Cancer Rep*. 2022;5(8).
42. Vinod R, Mahran R, Routila E, Leivo J, Pettersson K, Gidwani K. Nanoparticle-aided detection of colorectal cancer-associated glycoconjugates of extracellular vesicles in human serum. *Int J Mol Sci*. 2021;22(19).
43. Gidwani K, Huhtinen K, Kekki H, Van Vliet S, Hynninen J, Koivuviita N, et al. A nanoparticle-lectin immunoassay improves discrimination of serum CA125 from malignant and benign sources. *Clin Chem*. 2016;62(10):1390–400.
44. Hendrickson OD, Zherdev AV. Analytical Application of Lectins. *Crit Rev Anal Chem*. 2018;48(4):279–92.
45. Bojar D, Meche L, Meng G, Eng W, Smith DF, Cummings RD, et al. A Useful Guide to Lectin Binding: Machine-Learning Directed Annotation of 57 Unique Lectin Specificities. *ACS Chem Biol*. 2022;17(11):2993–3012.
46. Ward EM, Kizer ME, Imperiali B. Strategies and Tactics for the Development of Selective Glycan-Binding Proteins. *ACS Chem Biol*. 2021;16(10):1795–813.
47. Soukka T, Antonen K, Härmä H, Pelkkikangas AM, Huhtinen P, Lövgren T. Highly sensitive immunoassay of free prostate-specific antigen in serum using europium(III) nanoparticle label technology. *Clin Chim Acta*. 2003;212:45–85.
48. Boreham A, Brodwolf R, Walker K, Haag R, Alexiev U. Time-resolved fluorescence spectroscopy and fluorescence lifetime imaging microscopy for characterization of Dendritic polymer nanoparticles and applications in Nanomedicine. *Molecules*. 2017;22(1).
49. Härmä H, Soukka T, Lövgren T. Europium Nanoparticles and Time-resolved Fluorescence for Ultrasensitive Detection of Prostate-specific Antigen. *Clin Chem*. 2001;47(3):561–8.
50. Shen Y, Xu S, He D. A novel europium chelate coated nanosphere for time-resolved fluorescence immunoassay. *PLoS One*. 2015;10(6).
51. Yuan J, Matsumoto K, Kimura H. A New Tetradentate  $\beta$ -Diketone-Europium Chelate That Can Be Covalently Bound to Proteins for Time-Resolved Fluoroimmunoassay. *Anal Chem*. 1998;70(3):596–601.
52. Gidwani K, Kekki H, Terävä J, Soukka T, Sundfeldt K, Pettersson K. Nanoparticle-aided glycovariant assays to bridge biomarker performance and ctDNA results. *Mol Aspects Med*. 2020;72.
53. Kekki H, Peltola M, van Vliet S, Bangma C, van Kooyk Y, Pettersson K. Improved cancer specificity in PSA assay using Aleuria aurantia lectin coated Eu-nanoparticles for detection. *Clin Biochem*. 2017;50(1–2):54–61.
54. Ovarian cancer: recognition and initial management Clinical guideline [Internet]. [accessed on 20. 5. 2024]. Available from: [www.nice.org.uk/guidance/cg122](http://www.nice.org.uk/guidance/cg122)
55. Scholler N, Urban N. CA 125 in ovarian cancer. *Biomark Med*. 2007;1(4):513–23.
56. Saldova R, Struwe WB, Wynne K, Elia G, Duffy MJ, Rudd PM. Exploring the glycosylation of serum CA125. *Int J Mol Sci*. 2013;14(8):15636–54.
57. Wong NK, Easton RL, Panico M, Sutton-Smith M, Morrison JC, Lattanzio FA, et al. Characterization of the oligosaccharides associated with the human ovarian tumor marker CA125. *J Biol Chem*. 2003;278(31):28619–34.
58. Gidwani K, Huhtinen K, Kekki H, Van Vliet S, Hynninen J, Koivuviita N, et al. A nanoparticle-lectin immunoassay improves discrimination of serum CA125 from malignant and benign sources. *Clin Chem*. 2016;62(10):1390–400. »
59. Gidwani K, Nadeem N, Huhtinen K, Kekki H, Heinosaalo T, Hynninen J, et al. Europium Nanoparticle-Based Sialyl-Tn Monoclonal Antibody Discriminates Epithelial Ovarian Cancer-Associated CA125 from Benign Sources. *J Appl Lab Med*. 2019;4(3):299–310.
60. Jain S, Nadeem N, Ulfenborg B, Mäkelä M, Ruma SA, Terävä J, et al. Diagnostic potential of nanoparticle aided assays for MUC16 and MUC1 glycovariants in ovarian cancer. *Int J Cancer*. 2022;151(7):1175–84.

61. Salminen L, Nadeem N, Rolfsen AL, Dørum A, Laajala TD, Grønman S, et al. Exploratory Analysis of CA125-MGL and-STn Glycoforms in the Differential Diagnostics of Pelvic Masses. *J Appl Lab Med*. 2020;5(2):263–72.
62. Salminen L, Nadeem N, Jain S, Grønman S, Carpén O, Hietanen S, et al. A longitudinal analysis of CA125 glycoforms in the monitoring and follow up of high grade serous ovarian cancer. *Gynecol Oncol*. 2020;156(3):689–94.
63. Amiri M, Asadi Samani L, Kashi AH, Khadem N, Ziaee SAM, Mowla SJ. New Molecular Markers for Prostate Cancer Diagnosis. *UrolJ*. 2024;21(1):1–13.
64. Peracaula R, Tabarés G, Royle L, Harvey DJ, Dwek RA, Rudd PM, et al. Altered glycosylation pattern allows the distinction between prostate-specific antigen (PSA) from normal and tumor origins. *Glycobiology*. 2003;13(6):457–70.
65. Gabriele C, Prestagiacomo LE, Cuda G, Gaspari M. Mass spectrometry-based glycoproteomics and prostate cancer. *Int J Mol Sci*. 2021;22(10).
66. Song E, Mayampurath A, Yu CY, Tang H, Mehref Y. Glycoproteomics: Identifying the glycosylation of prostate specific antigen at normal and high isoelectric points by LC-MS/MS. *J Proteome Res*. 2014;13(12):5570–80.
67. Duffy MJ, Evoy D, McDermott EW. CA 15-3: Uses and limitation as a biomarker for breast cancer. *Clin Chim Acta*. 2010;411(23–24):1869–74.
68. Terävä J, Tiainen L, Lamminmäki U, Kellokumpu-Lehtinen PL, Pettersson K, Gidwani K. Lectin nanoparticle assays for detecting breast cancer-associated glycovariants of cancer antigen 15-3 (CA15-3) in human plasma. *PLoS One*. 2019;14(7).
69. Chandrasekar D, Guerrier C, Alisson-Silva F, Dhar C, Caval T, Schwarz F, et al. Warning signs from the crypt: Aberrant protein glycosylation marks opportunities for early colorectal cancer detection. *Clin Transl Gastroenterol*. 2023;14.
70. Fan J, Chen B, Luo Q, Li J, Huang Y, Zhu M, et al. Potential molecular biomarkers for the diagnosis and prognosis of bladder cancer. *Biomed Pharmacother*. 2024;173.
71. Ideo H, Kondo J, Nomura T, Nonomura N, Inoue M, Amano J. Study of glycosylation of prostate-specific antigen secreted by cancer tissue-originated spheroids reveals new candidates for prostate cancer detection. *Sci Rep*. 2020;10(1).
72. Tai CS, Lin YR, Teng TH, Lin PY, Tu SJ, Chou CH, et al. Haptoglobin expression correlates with tumor differentiation and five-year overall survival rate in hepatocellular carcinoma. *PLoS One*. 2017;12(2).